



- Host weight 1.4KG, Lightweight and facility
- HD large color screen Monitoring waveform dual interface display
- better fan quiet less than 30dB
- Integrated water tank Anti-backflow solution

Parameters

Type	CPAP			BIPAP		
Model (3.5 inch)	F-20C	F-20i	F-20A	F-25V	F-30S	F-30V
Model (5 inch)	R-20C	R-20i	R-20A	R-25V	R-30S	R-30V
Work Mode	CPAP	CPAP APAP iAPAP	CPAP, APAP	CPAP, S/T S, T, APCV VAF-ST VAF-APCV	CPAP S/T, S, T APCV, APAP iAPAP, AutoB	CPAP, S/T S, T, APCV VAF-ST VAF-APCV HF
Pressure Range	4-20	4-20	4-20	4-25	4-30	4-30
Boost Sensitivity	/	3 Levels	3 Levels	/	3 Levels	/
SPO ₂ Kit	Optional					
Heating tube	Optional					
Data management cloud platform	•					
Ramp	0-60 minutes					
Automatic switch	•					
Automatic leak compensation	•					
Automatic altitude adjustment	•					
Noise	<30 dB					
Humidifier	0-5 Level					
Size	243mm*181mm*110mm					
Weight	1370g					

ResFree
CPAP/BIPAP

- R-20C R-20i R-30A R-25V R-30S R-30V
- F-20C F-20i F-30A F-25V F-30S F-30V

Professional intelligent | comfort treatment


Intelligent constant temperature preheating in advance

Optional function of heating pipeline can ensure the output gas temperature and make the patient more comfortable during the treatment. Preheat in advance to avoid using the air too cold in winter.



Automatic light recognition

Under different light intensities, the screen brightness is automatically adjusted, which is clearly visible under strong light and away from the stimulation of bright lights at night.



24 hour service hotline // +86 18674382266 (WeChat & WhatsApp)

Email : sale@csbeyond.com URL : www.hnbyond.com

Add : Beyond Zone, Lijiacun Rd, Xueshi Street, Yuelu District, 410208, Changsha, China



Scan for video

